**911**  
**NETFULL**  
# 880020000KK

**KONG spa**

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## INDEX

<b>1 GENERAL INFORMATION</b>	<b>Page 3</b>
<b>2 TECHNICAL FEATURES</b>	<b>Page 4</b>
2.1 Terminology of parts and materials	Page 4
2.2 Accessories	Page 5
2.3 Dimensions	Page 6
2.4 Capacity	Page 7
<b>3 SPECIFIC INFORMATION</b>	<b>Page 8</b>
3.1 General	Page 8
3.2 Stretcher assembly	Page 8
3.3 Positioning and securing the patient	Page 10
3.4 Patient transport	Page 10
3.4.1 Connecting the handles	Page 11
3.4.2 Connecting the hanging kit	Page 11
<b>4 MAINTENANCE AND REPAIR</b>	<b>Page 12</b>
4.1 General	Page 12
4.2 Maintenance	Page 12
4.3 Repairs	Page 12
<b>5 STORAGE</b>	<b>Page 13</b>
<b>6 INSPECTIONS AND AUDITS</b>	<b>Page 14</b>
6.1 Inspections before and after use	Page 14
6.2 Audits	Page 14
<b>7 PRODUCT LIFE AND WARRANTY</b>	<b>Page 15</b>
7.1 Product life	Page 15
7.2 Warranty	Page 15
7.3 Legal obligations	Page 15

## 1 – GENERAL INFORMATION

Information provided by the manufacturer (the following information) must be read and understood by the user prior to use. This information relates to the description of the features, performance, assembly, disassembly, maintenance, care, disinfection, etc. of the device. Even if it contains some suggestions of use, it should not be considered a users manual in real situations.

### Warning:


- this device must only be used by persons who are physically able, informed and trained in its use and with specific experience in patient handling, or it may be used in training activities by persons under the direct supervision of trainers / supervisors that guarantee safety,
- all checks described in paragraph 6 should be carried out before and after use. If the user has the slightest doubt about the efficiency of the device, it must be replaced immediately,
- to reduce the risks of exposure / transmission of infectious diseases, the device must be cleaned and disinfected as described in paragraph 4,
- carefully follow the manufacturer's information, improper use of this device is dangerous,
- incorrect use of patient restraint systems can put the patient's life at risk,
- use in combination with devices and / or accessories other than those listed in paragraph 2.2 can be dangerous. You should always check the compatibility of devices using the manufacturer's information,
- improper use, deformations, drops, wear, chemical contamination, exposure to temperatures below -30°C or higher than +50°C for the components/devices textiles/plastics, and +100°C for metallic devices, are some examples of causes that can reduce, limit and end the life of the device,
- prior to any rescue operation, make sure the weight does not exceed the capacity defined in paragraph 2.4, and that it is uniformly distributed,
- avoid exposing the device to heat sources, and avoid contact with chemicals. Reduce unnecessary exposure to direct sunlight. At low temperatures and in the presence of moisture, ice can form that may reduce the flexibility of textile devices and increase the risk of cutting and abrasion,
- it is absolutely forbidden to modify and / or repair this device.

All our devices are tested and inspected piece by piece in accordance with Quality System procedures certified according to UNI EN ISO 9001. **Attention: laboratory tests, testing, information and the standards do not always manage to reproduce practice, so the results obtained in real conditions of use, in a natural environment, may sometimes differ to a considerable degree. The best information is the continued practice of use under the supervision of a competent person / expert / qualified person.**

## 2 – TECHNICAL FEATURES

### 2.1 – Terminology of parts and materials (fig. 1)

- a) Aluminum alloy frame, can be dismantled into two symmetrical parts,
- b) Polypropylene net,
- c) Polyurethane foam mattress,
- d) Polyester patient restraint straps,
- e) Aluminum alloy handles,
- f) Polyurethane padding for handles,
- g) Polyamide hanging kit with aluminum alloy connectors.

## 2.2 – Accessories

- a) Spinal boards “X-TRIM” (fig. 2) - 85510N000KK - 85511N000KK - 85512N000KK


Fig. 2 – Article 85510N000KK

Fig. 2 – Article 85511N000KK


Fig. 2 – Article 85512N000KK

b) Vacuum mattress "VACUUM" (fig. 3) - 8720040000KK

Fig. 3 – Article 8720040000KK

### 2.3 – Dimensions

- Length with handles: 330 cm
- Length without handles: 220 cm
- Width: 60 cm
- Weight: 8 kg

## 2.4 - Capacity

The stretcher "911 NET FULL" is tested to support the following uniformly distributed loads:

- 1500 kg when lifted by a hanging kit,
  - 450 kg when lifted and carried by the handles, to which one must apply the following minimum safety ratios, considered in relation to the type of use:
- |   |                          |
|---|--------------------------|
| - stretcher lifting and transport by means of the handles (fig. 4): | 1:3 = user load: 150 kg  |
| - lifting, by hanging kit, with rescue devices (fig. 5):            | 1:10 = user load: 150 kg |
| - lifting and transport, by hanging kit, with helicopters (fig. 5): | 1:14 = user load: 100 kg |

Fig. 4


Fig. 5

### 3 – SPECIFIC INFORMATION

#### 3.1 - General

The **880.020** stretcher “**911 NET FULL**” is a device suitable for recovering and transporting a patient, also while immobilizing them on the “X-TRIM” spinal board or on the “VACUUM” vacuum mattress.

The hanging kit permits stretcher lifting, lowering and transport using rescue devices and helicopter winches.

#### 3.2 – Stretcher assembly

Assemble the two parts of the structure as follows:

- a) Rotate the two safety rods (fig. 6)

Fig. 6

- b) Lift and rotate the locking rods (fig. 7),

Fig. 7

- c) Connect the two parts of the structure by joining the frame pins (fig. 8),

Fig. 8

- d) Reposition the two locking rods into their respective seats (fig. 9),

Fig. 9

- e) Reposition the two safety rods (figs. 10-11) and check the tightness of the connection,

Fig. 10

Fig. 11

- f) Fully screw the two sleeves into place (fig. 12).


Fig. 12

### 3.3 – Positioning and securing the patient

With the stretcher placed on a flat, stable surface:

- a) unclip the buckles and open the restraint straps,
- b) place the patient on the stretcher, clip the buckles on the restraint straps and tension the straps until they lock (figs. 13-16).

Fig. 13

Fig. 14

Fig. 15

Fig. 16

This procedure must also be applied in cases where the patient is immobilized using the “X-TRIM” spinal board or the “VACUUM” vacuum mattress.

### 3.4 – Patient transport

The stretcher can lift and transport the patient by means of:

- a) the lateral part of the frame,
- b) the handles,
- c) the hanging kit.

### 3.4.1 – Connecting the handles

- a) Connect the handles to the stretcher frame (fig. 17), respecting the connecting numbers,

Fig. 17

- b) Lock the connections in place by inserting the contoured pin into its hole (fig. 18).


Fig. 18

### 3.4.2 – Connecting the hanging kit

- a) Connect the hanging kit to the appropriate side holes of the stretcher using the eight connectors (fig. 19),


Fig. 19

- b) Fully screw the connector sleeves closed.

## **4 – MAINTENANCE AND REPAIR**

### **4.1 - General**

The “911 NET FULL” stretcher is made of materials with a high resistance to wear and tear. Despite this, the conditions of use call for maintenance work and, in particular cases, also repair.

### **4.2 - Maintenance**

Maintenance work that must be performed by the user are:

- a) cleaning: wash with warm fresh water (max. 40°C) after each use, possibly with the addition of a neutral detergent (e.g. Marseille soap). Rinse and allow to dry in the shade, away from direct heat,
- b) to disinfect, when deemed necessary: dilute bleach (sodium hypochlorite) in warm water with the concentration of 1% bleach and soak the product for one hour. Then rinse thoroughly with clean water and allow to dry in the shade, away from direct heat.

### **4.3 - Repairs**

Repairs must be carried out exclusively by the manufacturer and recorded on the inspection sheet.

The user is permitted to replace:

- a) the handles,
- b) the hanging kit,
- c) the patient restraint straps.

## **5 - STORAGE**

After cleaning, any disinfection, and drying, store the stretcher and its accessories in a place which is dry (relative humidity 40-90%), cool (temperature 5-40°C) and dark (avoid U.V. rays), chemically neutral (absolutely avoid saline and/or acidic environments), away from sharp edges, heat, humidity, corrosive substances and other potentially harmful conditions. **Do not store these products wet!**

## **6 – INSPECTIONS AND AUDITS**

### **6.1 – Inspections before and after use**

To ensure product efficiency and the safety of the patient and rescue team, the stretcher and its accessories must be inspected before and after every use. In particular, it is necessary to check that:

- a) the textile components are not cut or torn, especially in areas that come into contact with the buckles, connectors and structure,
- b) the stitching does not have loose or worn threads,
- c) the metallic parts have no deformations, cracks or wear,
- d) the connectors are operating correctly,
- e) the rivets are not worn or cracked,
- f) the mesh of the net is not broken.

### **6.2 - Audits**

Unless the provisions of law are more restrictive, the manufacturer recommends this device be submitted to Kong S.p.A. for an audit every three years from the date of first use. With inspections performed by a person trained by KONG S.p.A. at least annually, the frequency of audits recommended is five years from the date of first use.

## **7 – PRODUCT LIFE AND WARRANTY**

### **7.1 – Product life**

The product life is 15 years from the date of first use, provided that:

- a) the maintenance and storage are carried out as described in paragraphs 4 and 5 respectively,
- b) the inspections before and after use and the the audits do not indicate operating defects, deformations, wear, etc.,
- c) the product is used correctly.

Isolate devices that have not passed before use, after use, and periodic inspections.

### **7.2 - Warranty**

The manufacturer guarantees that the device conforms to the applicable standards at the time of production. The warranty for defects is limited to defects in raw material and workmanship, it does not include normal wear and tear, oxidization, damages caused by improper use and / or competitions, from improper maintenance, transportation, care or storage, etc.. The warranty is immediately void in the event that any changes are made to the device and / or it is tampered with. The validity corresponds to the legal warranty of the country in which the device was purchased, with effect from the date of sale by the manufacturer. Upon expiry of this period no claims can be made against the manufacturer. Any request for repair or replacement under warranty must be accompanied by a proof of purchase. If the defect is recognized, the manufacturer undertakes to repair or, at its discretion, replace or refund the device. In no case will the liability of the manufacturer extend beyond the invoice price of the device.

### **7.3 – Legal obligations**

Professional activities and leisure activities are often regulated by specific national laws which may impose restrictions and / or requirements on the use of PPE and to the provision of safety systems of which PPE are components. It is the obligation of the user to know and to apply these laws that could provide different limits to what is reported in this information.

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