**LECCO 2.0**

**870.04**

Directive 93/42/EEC


**KONG**

**KONG**

**EN**

[WWW.KONG.IT](http://WWW.KONG.IT)

Laden Sie die Übersetzung in Ihrer Sprache herunter  
- Download the translation in your language - Bájate  
la traducción en tu idioma - Télécharger la traduction  
dans votre langue - Scarica la traduzione nella tua  
lingua - 下载您语言版本的手册

# TABLE OF CONTENTS

<b>1 - SYMBOLS AND SUPPORT</b>	<b>4</b>
1.1 Symbols	4
1.2 Support	4
<b>2 - GENERAL INFORMATION</b>	<b>5</b>
<b>3 - TECHNICAL FEATURES</b>	<b>6</b>
3.1 Terminology and materials of the parts	6
3.2 Dimensions	8
3.3 Capacity	8
3.4 Accessories and spare parts	9
<b>4 - SPECIFIC INFORMATION</b>	<b>16</b>
4.1 Intended use	16
4.2 Preparing the scoop stretcher	17
4.3 Positioning and immobilising the patient	20
4.4 Transporting the patient	22
4.5 Disassembling the scoop stretcher	25
4.6 Removing the sheet from the metal structure	26
4.7 Fixing the sheet into the metal structure	26
<b>5 - MAINTENANCE AND REPAIR</b>	<b>30</b>
5.1 General	30
5.2 Maintenance	30
5.3 Repair	30
<b>6 - STORAGE</b>	<b>31</b>
<b>7 - CHECKS, INSPECTIONS AND SERVICING</b>	<b>32</b>
7.1 Checks	32
7.2 Inspections	32
7.3 Servicing	32
<b>8 - DURATION AND WARRANTY</b>	<b>33</b>
8.1 Duration	33
8.2 Disposal	33
8.3 Warranty	33
8.4 Law obligations	33
<b>9 - LABELLING AND SYMBOLS</b>	<b>34</b>
9.1 Labelling	34
9.2 Symbols	34
<b>10 - DOCUMENTS</b>	<b>35</b>
10.1 Recording inspections and servicing	35
10.2 Declaration of conformity (facsimile)	36
<b>11 - LIST OF APPLIED AND REFERENCE STANDARDS</b>	<b>37</b>
11.1 Applied standards	37
11.2 Standards used as reference	37

### 1.1 SYMBOLS

For the sake of comfort and clarity while reading this manual, below are the symbols used for handling important warnings for a proper and safe use of the device.

#### **REQUIREMENT FOR PROPER USE**

It identifies the presence of information for proper use of the device.

---

#### **INFORMATION REQUIREMENT**

It identifies the presence of useful and general information which guides the reader towards a conscious use of the device and/or performance of actions.

---

It shows that the product is made, designed and manufactured in compliance with the safety requirements (RES) laid down by the Medical Devices Directive 93/42/CEE (Class I medical devices, in compliance with the classification regulation 1 as shown in annex IX).

---

### 1.2 SUPPORT

For information please contact Kong Customer Support Service by:

- telephone 0039 0341 630506
- fax 0039 0341 641550
- email: [safetycare@kong.it](mailto:safetycare@kong.it)

or write to KONG S.p.A. – Via XXV Aprile, 4 – 23804 Monte Marenzo LC - ITALY.

To facilitate support operations, please always communicate or state the serial number (SN) indicated on the label fixed to the Medical Device.

# CHAPTER 2

## 2 GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the device. The information regards the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers tips on use, this information shall not be considered as a user manual under actual conditions of use.

### **WARNINGS AND LIMITATIONS OF USE:**

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, falls, wear, chemical contamination, exposing textile/plastic components/devices

- to temperature below  $-30^{\circ}\text{C}$  or above  $+50^{\circ}\text{C}$  and metal components/devices to temperatures exceeding  $100^{\circ}\text{C}$  are some examples of causes that can reduce, limit and end the life of the device,
- Prior to any rescue operation, be keen not to exceed the capacity indicated in paragraph 3.3,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or accessories different from the ones indicated in paragraph 3.4 can be hazardous,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase cutting and abrasion-related risks on textile and synthetic devices.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device in the natural environment may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/expert/qualified people.

**Warning: not suitable for use in an ATEX environment (Directive 94/9/EC)**

# CHAPTER 3

## TECHNICAL FEATURES

The scoop stretcher is provided, and can be kept away in the user's bag, with the sheet connected to the framework so as to facilitate preparation to use.  
The parts, illustrated in paragraph 3.1, can still be separated from each other.

### 3.1 TERMINOLOGY AND MATERIALS OF THE PARTS

A – Aluminium alloy heads,  
B – Aluminium alloy longitudinal beams,  
C1 and C2 – Aluminium alloy handles,  
D – Polyamide handle padding.


Fig. 1


Fig. 1

## 3.2 DIMENSIONS

Length without handles:	198 cm
Length with handles:	310 cm
Width:	46.5 cm
Overall weight (fig.1):	16.9 kg

## 3.3 CAPACITY

The table below defines the capacities of the “LECCO 2.0” stretcher, based on the method of use deemed most appropriate by the rescue operator, the outcome of the laboratory tests and recommended safety coefficients.

How to use	Tensile strength	Safety coefficient	Capacity
<b>Transportation with handles</b> (fig. 22, 23 and 24) <i>including model 871.50 HANDLES and 871.51 GRIPS (accessories)</i>	> 850 kg	1:4	200 kg
<b>Lifting and lowering with rescue devices</b> (fig. 2 and 3) <i>with 871.13 ORION 2.0 and 871.04 EVOLUTION suspensions (accessories)</i>	> 2100 kg	1:10	200 kg
<b>Lifting and lowering using helicopters</b> (fig. 2 and 3) <i>with the 871.13 ORION 2.0 and 871.04 EVOLUTION suspensions (accessories)</i>	> 2100 kg	1:14	150 kg
<b>Transportation on ropeway conveyor and overhead lines</b> (fig. 4) <i>with the 871.01 GEMINI suspension (accessory)</i>	> 2100 kg	1:10	200 kg
<b>Lowering and rescue on areas with steep inclinations</b> (fig. 5) <i>with the 871.20 VIRGO belt (accessory)</i>	> 850 kg	1:4	200 kg
<b>Lowering and rescue in space</b> (fig. 5) <i>with the 871.20 VIRGO belt (accessory)</i>	> 1250 kg	1:6	200 kg

Prior to any rescue operation, be keen not to exceed the indicated capacity

## 3.4 ACCESSORIES AND SPARE PARTS

The methods of connection and use of the accessories are defined in the relative user instructions.

### 3.4.1 Accessories

#### 871.13 ORION 2.0

Suspension with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for lifting and lowering using rescue devices and helicopters.


Fig. 2

#### 871.04 EVOLUTION

Two-piece suspension with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for lifting and lowering using rescue devices and helicopters.


Fig. 3

### 871.01 GEMINI

Suspension with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for transporting the scoop stretcher on ropeway conveyor and overhead lines.


Fig. 4

### 871.02 VIRGO

High-load-capacity belt, to be connected to the rope, for fastening when lowering and rescuing on areas with steep inclinations and in the space.


Fig. 5

### 871.03 AQUARIUS

Kit for fixing the scoop stretcher to the deck of the helicopters. Consisting of 4 length-wise adjustable webbings with oval links and quick-release stainless steel connectors.


Fig. 6

### 871.20 WILLY

Height-wise adjustable wheel with low-pressure tyre to be applied to the scoop stretcher in a barycentric position. It allows the transportation of the patient on flat grounds or on mule tracks.


Fig. 7

### 871.22 TWIN WILLY

Height-wise adjustable twinned wheel with low-pressure tyres to be applied to the scoop stretcher in a barycentric position. It allows the transportation of the patient on flat grounds or on mule tracks. Can be combined with the 871.20 WILLY single wheel.


Fig. 8

### 871.31 VISI

Lexan® face protection visor.


Fig. 9

### 871.50 HANDLES

Quick-coupling handles for transporting the scoop stretcher using the 871.20 WILLY and 871.22 TWIN WILLY wheels. Also suitable for shoulder-transportation carried by four rescue operators.


Fig. 10

### 871.51 GRIPS

Quick coupling grippers for moving the scoop stretcher in narrow spaces.




Fig. 11

**871.41 WRAP**

Backbone protection for transportation and sliding on snowy slopes.

**Fig. 12**

## 3.4.2 Spare parts

C1 and C2 – Handles,  
D– Paddings for handles,  
E – Sheet,  
H – Headrest cushion,  
M – Bag.


Fig. 13

# CHAPTER 4

## SPECIFIC INFORMATION

### 4.1 INTENDED USE

The **“LECCO 2.0”** scoop stretcher is a medical device, particularly suitable for mountain rescue operations, for rescuing and transporting a patient, even immobilised on **“X-TRIM”** spinal boards or on **“VACUUM”** mattresses.

Decisions concerning moving and immobilising the patient, as well as the duration, method to be used and combination with other devices, shall be taken and executed by expert and trained personnel only.

The 871.13 ORION 2.0 and 871.04 EVOLUTION suspensions, not part of standard supply, make the scoop stretcher utilisable with winch.

## 4.2 PREPARING THE STRETCHER

- Take the pre-assembled scoop stretcher from the bag and spread it with the sheet facing towards the ground,
- connect the central longitudinal beams (A1) to the longitudinal beams of the head element (A2) and subsequently to the longitudinal beams of the other head element.

To connect the longitudinal beams:

- keep the pin (N) pressed,
- insert the central longitudinal beam into the longitudinal beam of the head element,
- check whether the pin (N) is automatically inserted into the hole of the longitudinal beam of the head element.


Fig. 14

The resistance of the scoop stretcher is guaranteed by the correct assembly and stretching of the sheet belts. The pins (N) are provided solely for facilitating assembly.

c) Slightly tighten the transversal belts,


Fig. 15

d) slightly tighten the longitudinal belts being kept not to misalign the transversal belts,


Fig. 16

e) tighten the transversal belts definitely, insert the surplus part into the through element and position it between the sheet and the stretched belt,


Fig. 17

f) stretch the longitudinal belts definitely, insert the surplus part into the through element and position it between the sheet and the stretched belt.


Fig. 18

## 4.3 POSITIONING AND IMMOBILISING THE PATIENT

With the scoop stretcher positioned on a flat and stable surface:

a) open the sheet zip fully and release the buckles of the internal belts,


Fig. 19

b) lay the patient on the scoop stretcher, couple the buckles and stretch the belts until the patient is immobilised,

This procedure shall be applied even when patients are immobilised in “X-TRIM” spinal boards (except for the groin belt) or with “VACUUM” mattresses.


Fig. 20

c) close the sheet zip, couple the external belts and stretch them for optimal immobilisation of the patient.


Fig. 21

## 4.4 TRANSPORTING THE PATIENT

The scoop stretcher, with provided handles, is suitable for “shoulder” transporting the patient (figures 22, 23 and 24). Furthermore, it is used with the following accessories:

- a) 871.50 HANDLES and 871.51 GRIPS (figures 10 and 11),
- b) 871.00 ORION, 871.01 GEMINI and 871.04 EVOLUTION suspensions (figures 2, 3 and 4),
- c) 871.02 VIRGO belt (figure 5),

other methods of transportation can be applied, lift it and lower it, as indicated in the specific user instructions.

### 4.4.1 Connection with provided handles

- a) Holding the automatic peg (O) pulled, insert the provided handles into the relative seats of the head elements,
- b) release the automatic pin, check whether it is well-fitted in the slot of the handle and locks the connection.


Fig. 22

The provided handles can be exchanged between (can be attached to the metallic structure both at head and feet side) to allow the transport of the patient with **a suitable arrangement according to the rescuer opinion.**

Figures 23, 24, and 25 are **only hint of handles configurations** to maintain the “head above feet” configuration while transporting the patient on different sloping grounds.


Fig. 23 - Transport on flat or slightly sloping grounds.


Fig. 24 - Downhill transport on steep and vertical slopes.


Fig. 25 – Uphill transpor.

## 4.5 DISASSEMBLING THE STRETCHER

With the stretcher upside-down (sheet facing upwards):

a) loosen the transversal and longitudinal belts enough to be able to disconnect the longitudinal beams,


Fig. 26

b) press the buttons (N) and disconnect the longitudinal beams of the head elements (A2) from the central longitudinal beams (A1),


Fig. 27

c) fold the stretcher into three parts and put it away in the bag.

## 4.6 REMOVING THE SHEET FROM THE METAL STRUCTURE

Whenever it is required to remove the sheet from the metal structure (e.g. for disinfection, substitution, etc.):

- a) turn upside down the stretcher (with the sheet above);
- b) loosen the transversal and longitudinal belts (fig. 26) and remove them from the through elements;
- c) remove the through holes from the millings of the metal structure.

The through elements should not be overlapped so as to be able to remove the through holes from the millings.

## 4.7 FITTING THE SHEET INTO THE METAL STRUCTURE

- a) connect the central longitudinal beams (A1) to the longitudinal beams of the head element (A2) and subsequently to the longitudinal beams of the other head element, matching the marked identification numbers (1-1, 2-2, etc);


Fig. 28

b) spread the sheet with the bottom facing upwards,

The through elements are positioned in an alternating fashion so as to provide the sheet-bearing function even in case of accidental detachment.


Fig. 29

c) overturn the metal structure and place it on the sheet so that the position of the millings perfectly matches the transversal belts,


Fig. 30

d) insert the through elements of the transversal belts, from the external towards the internal, in the relative millings of the metal structure.

Ensure that the through elements do not overlap, so as to be able to insert them.


Fig. 31

e) Insert the transversal belts into the through elements and tighten them slightly.


Fig. 32

f) Pass the longitudinal belts above the crosspieces of the head elements, insert them into the through elements and tighten them slightly, being keen not to misalign the transversal belts.

Proceed as indicated:

- from point e) of paragraph 4.2 to prepare the stretcher for use,
- in paragraph 4.5 to put the stretcher away into the bag.


Fig. 33

# CHAPTER 5

## MAINTENANCE AND REPAIR

### 5.1 GENERAL

The “**LECCO 2.0**” stretcher is made of material highly resistant to wear and external agents. Despite this, the use conditions require maintenance operations and even repair in particular conditions.

Maintenance and repair operations must be recorded in a special register an example of which is indicated in chapter 10.

### 5.2 MAINTENANCE

Below are maintenance operations to be carried out by the user:

- a) cleaning: after every use, wash with warm (max. 40° C) potable water, possibly adding a neutral detergent (e.g. Marseille soap). Rinse and leave to dry in a shaded place, away from direct sources of heat,
- b) disinfection, when deemed necessary: soak the product in warm water in which sodium hypochlorite has been diluted at 1% for one hour, then rinse thoroughly using potable water and leave to dry in a shaded place, away from direct sources of heat.

### 5.3 REPAIR

Repair shall be strictly carried out by the manufacturer.

The user is allowed to replace the parts mentioned in paragraph 3.4.2 only, with new and genuine parts.

# CHAPTER 6

## STORAGE

### 6 STORAGE

After cleaning, any disinfection and drying operations, keep the product and its accessories in a place that is dry (relative humidity 40-90%), cool (temperature 5-40° C) and safe (avoid UV radiations), chemically neutral (absolutely avoid salty and/or acid environments), away from sharp edges, sources of heat, humidity, corrosive substances or other possible jeopardising conditions.

**Do not store this device wet!**

# CHAPTER 7

## CHECKS INSPECTIONS AND SERVICING

### 7.1 CHECKS

Before and after every use, check the device and ensure that:

- it is appropriate for the use it was intended for,
- it has not been subjected to mechanical deformations and it neither has cracks nor shows signs of wear,
- the textile parts do not have cuts, burns, chemical products residues, excessive hair, wear, in particular check areas at contact with the metal components,
- markings, including labels, are readable.

### 7.2 INSPECTIONS

The device must be inspected at least once a year, starting from the date of first use, by expert personnel trained and approved by KONG S.p.A. The date of the first use and the outcome of the inspection must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

### 7.3 REVISIONS

The device must be serviced by KONG S.p.A., or by expert personnel trained and approved by KONG S.p.A. if:

- any malfunctions are detected,
- the result of the pre and post use check or of the inspections is negative.

The outcome of the servicing must be recorded in the Inspections and Servicing Register, regarding which an example is indicated in chapter 10.

The device subject of servicing has a one-year warranty from the date of servicing.

# CHAPTER 8

## DURATION AND WARRANTY

### 8.1 DURATION

Carefully read the “Warnings and limitations of use” paragraph in chapter 2. The duration of this device depends on the outcome of the at least one-yearly inspection (paragraph 7.2) and possible servicing.

**Important: any devices that have not passed the pre-use, post-use and periodic inspection must be eliminated and made unusable.**

### 8.2 DISPOSAL

Follow the rules in force in the country of use and the waste disposal procedures of the reference hospital facility for proper disposal.

### 8.3 WARRANTY

The manufacturer guarantees the compliance of the device with the regulations in force at the time of production. Warranty covering defects is limited to defects regarding raw materials and manufacturing defects. The warranty does not cover normal wear, oxidation, damage caused by non-compliant use and/or in competition, by improper maintenance, transportation, preservation or storage. Warranty will be immediately deemed null and void should the device be modified or tampered with. Validity corresponds to the legal validity of the country where the device was sold, as from the date of sale, by the manufacturer. No claims shall be raised against the manufacturer once the aforementioned period has expired. Any repair or replacement request under warranty must be accompanied by a proof of purchase. Should the defect be recognised, the manufacturer undertakes to repair or, replace or refund the device at his discretion. The manufacturer shall not be held liable beyond the invoice price of the device in any case whatsoever.

### 8.4 LAW OBLIGATIONS


Professional and free time activities are often regulated by special national laws which can impose limits and/or obligations concerning the use of these devices. The user has the obligation to know and apply such laws which could provide for limits different from those indicated in this information.

# CHAPTER 9

## LABELLING AND SYMBOLS

### 9.1 LABELLING

Labelling example.

Any changes to the positioning of the symbols does not change the content thereof.

### 9.2 SYMBOLS

Manufacturer's identification

Product identification code

Unique serial number

LLLLLL : production batch  
YY : year of manufacture  
XXXX : progressive number

Consult the user manual

Warning: safety-related information, refer to the user manual

Compliant with Council Directive 93/42/EEC et seq. Class I MD

# CHAPTER 10

## DOCUMENTS

### 10.1 INSPECTIONS AND SERVICING REGISTER

<b>STRETCHER – Class I Medical Device</b>				
<b>REF.</b>	870.04	<b>FORM</b>	LECCO 2.0	
<b>SN</b>		<b>MD registration n°</b>	1776887	
<b>Date of first use</b>				
<b>INSPECTION AND SERVICING</b>				
<b>DATE</b>	<b>I/S</b>	<b>Description</b>	<b>Outcome</b>	<b>Manager</b>

## 10.2 DECLARATION OF CONFORMITY (FACSIMILE)

This device comes along with its EC Declaration of Conformity whose drafted original copy is signed. In case of loss or should the packaging be damaged, you may request replacement at: [safetycare@kong.it](mailto:safetycare@kong.it) submitting the unique serial number [SN] indicated on the label of the device.

**KONG s.p.a.**

via XXV Aprile, 4 | 23804 Monte Marengo (LC) - ITALY | tel +39 0341.630506 | [info@kong.it](mailto:info@kong.it)

certified UNI EN ISO 9001

### DECLARATION OF CONFORMITY

**Manufacturer:**

Company name : **KONG S.p.A.**  
Registered office: **Via XXV Aprile, 4 - 23804 Monte Marengo (LC)**  
Operational offices: **Via XXV Aprile, 4 - 23804 Monte Marengo (LC)**  
VAT Code: **IT 00703180166**

**hereby declares that the:**

**Medical device:** **LECCO 2.0**  
**Class:** **I according to attachment IX regulation 1**  
**REF:** **870.04**  
**SN (serial number):** \_\_\_\_\_  
**Date of manufacture:** \_\_\_\_\_  
registered in the Ministry of Health data bank under number **1776887**,

**complies with the essential requirements laid down by the Italian Legislative Decree n° 46/97 et seq. implementing Medical devices**

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Revision form 0 dated 01/08/2016

# CHAPTER 11

## LIST OF APPLIED AND REFERENCE STANDARDS

### 11.1 APPLIED STANDARDS

- EN 1041:2008  
Information supplied by the manufacturer of medical devices
- EN ISO 14971:2012  
Medical devices - Application of risk management to medical devices
- EN ISO 15223-1:2016  
Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
- EN 62366:2008  
Medical devices - Application of usability engineering to medical devices

### 11.2 STANDARDS USED AS REFERENCE


- EN 1865-1:2010+A1:2015  
Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
- EN ISO 13485:2016  
Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 13718-1:2008  
Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances

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